


# Seal Label Package - Set Hidrocor Toric



## Key list

- |  |   |
|--|---|
| <ul style="list-style-type: none"> <li>1 - Product and color</li> <li>2 - Catalogue number</li> <li>3 - Product 1<br/>PWR (Power), BC (Base Curve) and DIA (Diameter)</li> <li>4 - Product 1<br/>Serial Number</li> <li>5 - Batch code</li> <li>6 - Use-by date</li> </ul> | <ul style="list-style-type: none"> <li>7 - Product 2<br/>Serial Number</li> <li>8 - Product 2<br/>PWR (Power), BC (Base Curve) and DIA (Diameter)</li> <li>9 - EAN Bar</li> <li>10 - Unique device identifier:<br/>(01) EAN, (21) ID, (17) Expiration date and (10) Lot</li> <li>11 - Date of manufacture</li> <li>12 - custom-made device</li> </ul> |
|--|---|

## Importer and Exporter Label

1

2


**Solótica Indústria e Comércio Ltda.**  
CNPJ: 61.406.203/0001-77  
Rua Domingos Fasolari, 110 - Casa Verde  
São Paulo - SP. Tel: +55 11 3959-3500  
R.T.: Patrícia M.M. Duque CRF 27537  
Registro Anvisa: 10231649004

2024/07 Rev.00

### Key list

- 1 - Icons - Importer and Distributor
- 2 - Name and address of the importing and/or distributing entity

# Hidrocor Toric - Label



## Key list

- |   |  |
|---|--|
| 1 - Color   | 9 - Use-by date  |
| 2 - Base Curve and Diameter   | 10 - Content   |
| 3 - Batch code  | 11 - Power   |
| 4 - Icons - Consult instructions for use or consult electronic instructions for use / Do not resterilize/ Do not use if package is damaged and consult instructions for use / Single patient multiple use | 12 - Catalogue number  |
| 5 - Sterility Indication  | 13 - Authorized representative in the European Community/European Union                |
| 6 - EAN Bar   | 14 - Revision number and date  |
| 7 - Sterilized using steam heat Single sterile barrier system   | 15 - custom-made device  |
| 8 - Date of manufacture   | 16 - Manufacturer  |
|   | 17 - CE identification   |
|   | 18 - Unique device identifier:<br>(01) EAN, (21) ID, (17) Expiration date and (10) Lot |
|   | 19 - ID product  |