

Importer and Exporter Label

1

2




Solótica Indústria e Comércio Ltda.
CNPJ: 61.406.203/0001-77
Rua Domingos Fasolari, 110 - Casa Verde
São Paulo - SP. Tel: +55 11 3959-3500
R.T.: Patrícia M.M. Duque CRF 27537
Registro Anvisa: 10231649004

2024/07 Rev.00

Key list

- 1 - Icons - Importer and Distributor
- 2 - Name and address of the importing and/or distributing entity

Seal Label Package - Set Natural Colors




Key list

- 1 - Product and color
- 2 - Catalogue number
- 3 - Product 1
PWR (Power), BC (Base Curve) and DIA (Diameter)
- 4 - Product 1
Serial Number
- 5 - Batch code
- 6 - Use-by date

- 7 - Product 2
Serial Number
- 8 - Product 2
PWR (Power), BC (Base Curve) and DIA (Diameter)
- 9 - EAN Bar
- 10 - Unique device identifier:
(01) EAN, (21) ID, (17) Expiration date and (10) Lot
- 11 - Date of manufacture

Natural Colors - Label



Key list

- | | |
|--|---|
| <ul style="list-style-type: none"> 1 - Color 2 - Base Curve and Diameter 3 - Batch code 4 - Icons - Consult instructions for use or consult electronic instructions for use / Do not resterilize/ Do not use if package is damaged and consult instructions for use / Single patient multiple use 5 - Sterility Indication 6 - EAN Bar 7 - Sterilized using steam heat Single sterile barrier system 8 - Date of manufacture | <ul style="list-style-type: none"> 9 - Use-by date 10 - Content 11 - Power 12 - Catalogue number 13 - Authorized representative in the European Community/European Union 14 - Revision number and date 15 - Manufacturer 16 - CE identification 17 - Unique device identifier:
(01) EAN, (21) ID, (17) Expiration date and (10) Lot 18 - ID product |
|--|---|