

Seal Label Package - Unit Hidrocor Toric



Key list

- 1 - Product and color
- 2 - Catalogue number
- 3 - Product 1
PWR (Power), BC (Base Curve) and DIA (Diameter)
- 4 - Product 1
Serial Number
- 5 - Batch code
- 6 - Use-by date

- 7 - EAN Bar
- 8 - Unique device identifier:
(01) EAN, (21) ID, (17) Expiration date and (10) Lot
- 9 - Date of manufacture
- 10 - custom-made device

Importer and Exporter Label

1

2


**Solótica Indústria e Comércio Ltda.**
CNPJ: 61.406.203/0001-77
Rua Domingos Fasolari, 110 - Casa Verde
São Paulo - SP. Tel: +55 11 3959-3500
R.T.: Patrícia M.M. Duque CRF 27537
Registro Anvisa: 10231649004

2024/07 Rev.00

Key list

- 1 - Icons - Importer and Distributor
- 2 - Name and address of the importing and/or distributing entity

Hidrocor Toric - Label



Key list

- | | |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------|
| 1 - Color | 9 - Use-by date |
| 2 - Base Curve and Diameter | 10 - Content |
| 3 - Batch code | 11 - Power |
| 4 - Icons - Consult instructions for use or consult electronic instructions for use / Do not resterilize/ Do not use if package is damaged and consult instructions for use / Single patient multiple use | 12 - Catalogue number |
| 5 - Sterility Indication | 13 - Authorized representative in the European Community/European Union |
| 6 - EAN Bar | 14 - Revision number and date |
| 7 - Sterilized using steam heat
Single sterile barrier system | 15 - custom-made device |
| 8 - Date of manufacture | 16 - Manufacturer |
| | 17 - CE identification |
| | 18 - Unique device identifier:
(01) EAN, (21) ID, (17) Expiration date and (10) Lot |
| | 19 - ID product |